



[info@transomic.com](mailto:info@transomic.com)  
601 Genome Way, Suite 2021  
Huntsville, AL 35806  
866-833-0712

### **Safety and handling of lentiviral particles**

Recombinant lentivirus is considered a Biosafety Level 2 organism by the National Institutes of Health and the Center for Disease Control and Prevention. However, local health and safety regulations should be determined for each institution.

For more information on Biosafety Level 2 agents and practices, download Biosafety in Microbiological and Biomedical Laboratories (BMBL), Fifth Edition (Revised December 2009) published by the U.S. Department of Health and Human Services Centers for Disease Control and Prevention and NIH. The publication can be found here: <http://www.cdc.gov/biosafety/publications/bmb15/>.

If additional measures are needed, review biosafety guidance documents such as the NIH's "Biosafety Considerations for Research with Lentiviral Vectors" which refers to "enhanced BL2 containment". More information can be found through the NIH Office of Biotechnology Activities web site ([https://osp.od.nih.gov/wp-content/uploads/2014/01/Lenti\\_Containment\\_Guidance\\_0.pdf](https://osp.od.nih.gov/wp-content/uploads/2014/01/Lenti_Containment_Guidance_0.pdf))

### **Summary of Biosafety Level 2 Practices**

The following is meant to be a summary of Biosafety Level 2 practices and should not be considered comprehensive. A full account of required practices should be determined for each institute and/or department.

#### **Standard microbiological practices**

- Limit access to work area
- Post biohazard warning signs
- Minimize production of aerosols
- Decontaminate potentially infectious wastes before disposal
- Use precautions with sharps (e.g., syringes, blades)
- Review biosafety manual defining any needed waste decontamination or medical surveillance policies

#### **Safety equipment**

- Biological Safety Cabinet, preferably a Class II BSC/laminar flow hood (with a HEPA microfilter) used for all manipulations of agents that cause splashes or aerosols of infectious materials; exhaust air is not recirculated
- Protective personal equipment includes: protective laboratory coats, gloves, face protection if needed

#### **Facilities**

- Autoclave available for waste decontamination
- Chemical disinfectants available for spills